

PHS-3C RECHARGEABLE PH METER



1. Use of instruments

PHS type digital display acidity meter is a high-precision instrument used to measure PH value in laboratory. With dual functions of PH and MV, and standard e201-c composite electrode, it has quick reaction and good stability. The instrument is applicable to industrial and mining enterprises, scientific research, agriculture, medicine, teaching, environmental protection and other units, used for sampling and determination of pH value of aqueous solution and measurement of electrode potential (mV value).

2. Technical performance and complete set of instruments

Main technical parameters of the instrument:

Instrument model

PHS – 25

PHS – 2 – c

PHS – 3 – c

PH measurement range

The pH of 0.00 to 14.00

The pH of 0.00 to 14.00

The pH of 0.00 to 14.00

Accuracy of pH measurement

& plusmn; 0.05 pH

& plusmn; 0.02 pH

& plusmn; 0.01 pH

The resolution of the

0.01 pH

0.01 pH

0.01 pH

MV range of measurement

0 ~ & plusmn; 1999 mv