

WSB-2Y DESKTOP FLUORESCENT WHITENESS METER



Main technical parameters

Measurement range: 0 ~ 199

Content of measurement: blue ray whiteness r457 fluorescent whiteness F

Light source: LED light source

Lighting conditions: comply with GB/T3978 regulation 45/0

Display: LCD display

Minimum: 0.1

Measurement repeatability: 0.2

Error of indication: $\leq \pm 1$

Using the environment temperature (5 °C ~ 35 °C); Humidity ($\leq 85\%$ RH)

Transportation and storage environment: temperature (20 °C to 50 °C); Humidity ($\leq 90\%$ RH)