

KD-BC Tissue Conserving Table



FEATURES:

- KD-BC conserve any-type tissue embedding cassette, embedding box, embedding mold and tissue which was processed; high-capacity easy for storage, this instrument can be operated separately, also used for tissue spreading.

Technical specifications

Setting temperature	ambient ~90°C
Protection requirement of time delay	
Volume	4.5 liter
Voltage	AC 220V 50 Hz AC 110V 60Hz
Power	180W
Dimension	620×380×425 mm (W×D×H) (KD-BC)
Weight	15kg