

KD-BL Cooling Plate



Features

- KD-BL Cryo-plate adopted with new-type inverter compressor, which makes tissue paraffin cooled quickly Teflon plating ensure easy cleaning cooling temperature can be adjusted freely, temperature is settled properly to avoid tissue wax from not cooling enough or frost splitting, KD-BL cryo-plate can be use separately.

Technical specification

Area of specification	320x300mm
Voltage	AC220V 50Hz AC110V 60Hz
Power	300W
Dimension	590×345×385 mm (W×D×H), (KD-BL)
Weight	27kg