


KD-BLIII Cooling Plate



Features

KD-BLIII Cryo-plate adopted with new-type inverter compressor, which makes tissue paraffin cooled quickly Teflon plating ensure easy cleaning cooling temperature can be adjusted freely, the temperature is settled properly to avoid tissue wax from not cooling enough or frost splitting, KD-BL cryo-plate can be used separately.

Technical specifications

| | |
|--------------------------------------|------------------------------------|
| Setting temperature | RT~-20°C |
| Protection requirement of time delay | |
| Area of cryo-plate | 30x322mm |
| Voltage | AC220V 50Hz AC110V 60Hz |
| Power | 300W |
| Dimension | 685×380×430 mm (W×D×H) (KD-BL III) |
| Weight | 26kg |