


KD-BLIV Cooling System



Technical Specifications

Temperature Ranges of Cryo Module	0 to -20°C
Temperature setting of Cryo Module	the optimal working temperature is approximately -10°C, and the temperature is under delay protection
Dimensions of Cryo Module	315×380 mm
Working Voltage	AC 220V±10% 50Hz (standard model); AC110V±10% 60Hz
Power	300W (KD-BLIV)
Dimensions	710×350×390mm (W×D×H) (KD-BLIV)
Net weight	25kg (KD-BLIV)