

KD-BMR Paraffin Dispenser



FEATURES:

- Utilizes full computer automated control and intelligent digital temperature control
- Utilizes new type of heating elements that enable fast,, energy efficient and reliable heating
- Temperature sensors are made in USA, precise and reliable
- Automatically saves the latest operation settings
- Memory and automatic restoration function enables the preset temperature to be automatically saved after device operation
- LED displays the heating status

Major Technical Specification

Temperature Range	Continuously adjustable within 0 – 100°C
Temperature control precision	±1°C
Capacity	10000 ml
Working Voltage	AC 220V±10% 50Hz (standard model); AC110V±10% 60Hz
Power	1200W
Dimensions	355×410×540mm (W×D×H)
Net weight	16kg