

# KD-RS Fully Automated Tissue Stainer



## Features

- Integrated imported high-quality LCD touch-screen with intuitive software offers clear display, prompt response, high efficiency, and simple operation
- Unique ergonomic design of the staining protocol meets different users' needs
- Low noise, wear-resistant flexible transmission system using imported high-quality materials and elements to ensure reliable performance
- Intelligent automatic water influx/efflux/drainage system ensures the cleanness of water used in each step and improves water efficiency
- Special staining basket allows prompt staining of dozens of slides, which is a safe, reliable and green process
- Slides are rinsed and cleaned with circulating water, and it has the drying function.

## Technical Specification

- Number of Processing Cups: 14 (12 cups for reagents, the one at the 8th station for washout, and the one at the 1st station for drying)
- Length of Processing Time for each station: Adjustable within 0 min 0 s – 59 min 59 s for each station

Capacity of Single Cup	1500ml
Number of slides to be processed at the same time	72pcs
Working Voltage	AC 220V $\pm$ 10% 50Hz (standard model); AC110V $\pm$ 10% 60Hz
Power Requirements	500W
Dimensions	1180×420×470 mm (W×D×H)
Net weight	70kg