

KD-THII Floatation Workstation



FEATURES:

- Pure-green digital screen for real-time monitoring of heating temperature, direct and clear, easy to operate; all parameters including preset temperatures, working temperatures, and working status are real-time displayed
- Temperature is automatically program-controlled by single-chip microprocessors made in USA.
- This device features multiple functions and its easy setting operation can meet the needs of different users
- DC low-voltage illuminating system and removable transparent heating dish: easy operation and convenient observation
- New heating mechanism using a new-type high-thermal-conductivity heating element provides even and quick heating
- PID-controlled
- Triple temperature controls
- long lifespan, safe, reliable and energy-saving
- Temperature is precisely measured by temperature-sensing integrated blocks made in USA, and all settings are automatically stored in the system
- Special temperature-control and monitoring system characterized by strong resistance to abrasion and corrosion to ensure temperature precise
- Automatic memory and restoration functions: After startup, all preset temperature data are automatically stored in the system
- Temperatures for floating and drying slides are individually controlled
- The unit has the sounds system, it can put TF card and has UBS socket, also can use the headphones