



WXG-4_LC & WXG-4_ELC Polarimeter (Sodium & LED)

It is simple and useful Manual Polarimeter. Our goal is to be the first in World. Now we are doing more to do our best for you.

Features:

1. Easy to adjust scale with Streamline Body
2. Improve measurement error for optimized mould structure and precision scale
3. Delay operation time

Specifications

Range: $\pm 180^\circ$

Scale Division: 1°

Dial Venier Reading: $\pm 0.05^\circ$

Magnifier: 3X

Light Source: Sodium Light 589.44nm or LED Light

Test Tube: Round Cuvette with bubble trap and cap

Power Supply: 220V/50Hz or 110V/60Hz

Working Current: 1.3A

Discharging Power: 20W

Stabilization Time: 10minutes

LED



Sodium Light

