K37-24 ID Card Incubator



Product Information

The K37-24 ID-Incubator holds 2 independent incubation zones with a capacity for 12 Gel cards, and user can set temperature and incubator time.

It can be suitable many types ID-CARDS from different manufacturers.

1. Capacity: 24 Gel cards

2. Temperature: 37°C

3. Timer: 15min (Variable)

4. Power: AC220V or AC120V 50-60Hz, 300W

5. Dimensions: DxWxH=380x280x140 mm

Features

- 1. Incubator for 24-Gel Cards in racks at 37°C
- 2. 2 independent incubation zones with a capacity for 12 Gel cards
- 3. Digital interface screen for control and visualization of remaining incubation time and real temperature for

each zone

- 4. Fixed pre-set temperature of 37 degree
- 5. Pre-set incubation time of 15 minutes, which can be adjusted by the user
- 6. Audible programmable alarm to notify end of incubation period

Technical Data

Capacity	24 Gel cards
Temperature	37°C
Temperature uniformity	≤0.5°C
Temperature stability	≤0.2°C
Temperature accuracy	<0.5°C
Timer	15min (Variable)
	AC220V 3A 300W
Power	
Dimensions	380x220x140mm (DxWxH)