

MTV-1 Vortex



Product Information

The Mini vortex is the ideal vortex mixer for every lab with its compact size, ease of use, and its economic price point. The 2,500rpm fixed speed ensures there is no sacrificing any mixing power. It requires little working space and is easy-to-use.

Features

1. Touch activated vortex mixer
2. Stable compact design
3. 2,500rpm fixed vortexing speed
4. Powerful vortexing for tubes up to 15ml
5. MTV-1 is with a power adapter and a battery case
6. Choice of different colored housing

Technical Data

Model	MTV-1
Working style	Continuous/Short-Spin
Shaking movement	Orbital, 3mm
Speed Range	800~2500rpm
Dimensions(mm)	74 x 91 x 84
Voltage	DC6V, (4pcs 1.5V double A battery)
Weight	0.35Kg