## MX100-4A

**Product Information** 



The MX100-4A Is designed for shaking ELISA plate (96 well plates) TCPS and 1.5/2.0ml centrifuge tubes.

## Features

- 1. Variable speed from 200rpm to 1,500rpm
- 2. DC brushless motor for gentle and powerful vibration
- 3. Continuous or timing operation
- 4. Orbital of 3mm for efficient mixing in small sample vessels
- 5. Simple spring holders for easy and quick plate installation

## **Technical Data**

Model	MX100-4A
Mixing Orbit	3mm
Speed Setting Resolution	10rpm increment
Speed	200~1500 rpm @Standard 200~1000 rpm @Deep well plate or tube holder
Timer	1min~99h59min
Dimensions(mm)	AC220V or AC120V
Dimensions(mm)	290 x270 x 120
Net Weight	6.5kg

## Accessories

NO	Code	Description
01	AS-04040-00	MX100-4A Microplate shaker
02	AS-02041-01	20×1.5ml tube holder (apolegamy)
03	AS-02041-02	20×0.5m tube holder (apolegamy)