

Mini-15K Mini Centrifuge



Features

The Mini-15K is a practical, compact and reliable laboratory assistant which the highest speed reach 14,500rpm. Humanized design, simple and compact appearance, special processing of low noise design creates a quiet, comfortable experimental environment for you. Sensitive and efficient speed control, accelerated to the highest speed or reduce from the highest speed only 15sec.

Technical Data

Model	Mini-15K
Power requirement	100W
Relative centrifugal force (RCF)	14,000 x g
Max. rotational speed	14,500rpm
Timer	20sec~99min

Acceleration time to max. speed	15s
Braking time from max. speed	15s
Dimensions (DxWxH)	270x230x155mm

Accessories

Code	Description
AS-08041-01	Adapters for 0.5ml tubes
AS-08041-02	Adapters for 0.2ml tubes