

# Mini-6KC Mini Centrifuge



## Product Information

1. Mini-6KC Mini-Centrifuge is compact and cost effective, allowing any workstation to be equipped with a “personal” mini centrifuge.
2. It includes 2 quick release interchangeable rotors, two speed options (4,000rpm & 6000rpm).
3. The electronic brake can reach quick deceleration to reduce handling time.

## Features

1. Composite shell, shatter-resistant
2. The electronic braking system for rapid sample processing
3. Quick release rotor tub, allows for easy rotor exchange
4. Mini-6KC / Mini-6KS / Mini-6KT  
Two speed options: 4,000 & 6,000rpm
5. include 2 quick release interchangeable rotors  
–8×1.5/2.0ml tube rotor, –4 place PCR Strip rotor
6. 3-in-1 rotor design, high capacity and time-saving  
–6×1.5/2.0ml tubes, –6×0.5ml tubes, –2-8×0.2ml PCR strip tubes

## Technical data

Model	Mini-6KC
Capacity	8×1.5/2.0 ml tube rotor & 4 place PCR Strip rotor
Two Speed options	4,000rpm & 6,000rpm
Maximum RCF	1,200 & 2,000 x g
Dimensions	135 x 147 x 107mm (D x W x H)
Weight	0.7kg
Electrical	AC100~240V, 50/60Hz

## Accessories

Code	Description
AS-08020-00	Mini-6KC AC100V~240V 50~60Hz 1A
AS-08021-01	8×1.5/2.0ml rotor
AS-08021-02	4x8x0.2ml strip rotor
AS-08061-05	Adapters for 0.5/0.6ml tubes
AS-08061-06	Adapters for 0.2ml tubes