

# Mini-6KS Mini Centrifuge



## Product Information

1. Mini-6KS Mini-Centrifuge is compact and cost-effective, allowing any workstation to be equipped with a “personal” mini centrifuge.
2. It is supplied with one rotor suitable for 6-1.5/2.0ml tubes, 6-0.5ml tubes and 2-8×0.2ml PCR tubes, and two-speed options (4,000 & 6,000rpm).
3. The electronic braking can reach quick deceleration to reduce handling time.

## Feature:

1. Mini-6KS Mini-Centrifuge is compact and cost effective, allowing any workstation to be equipped with a “personal” mini centrifuge.
2. It is supplied with one rotor suitable for 6-1.5/2.0ml tubes, 6-0.5ml tubes and 2-8×0.2ml PCR tubes, and two speed options (4,000 & 6,000rpm).
3. The electronic braking can reach quick deceleration to reduce handling time.
4. Composite shell, shatter-resistant
5. The electronic braking system for rapid sample processing
6. Quick release rotor tub, allows for easy rotor exchange
7. Two speed options: 4,000 & 6,000rpm
8. include 2 quick release interchangeable rotors-8×1.5/2.0ml tube rotor, –4 place PCR Strip rotor

## Technical Data

Model	Mini-6KS
Capacity	6×1.5/2.0ml, 6×0.5ml tubes and 2-8×0.2ml PCR Strip rotor
Two Speed options	4,000rpm & 6,000rpm
Maximum RCF	1,200 & 2,000 x g
Dimensions	135 x 147 x 120mm (D x W x H)
Weight	0.7kg
Electrical	AC100~240V, 50/60Hz

## Accessories:

Code	Description
AS-08031-01	6×1.5/2.0ml, 6×0.5ml tubes and 2-8×0.2ml PCR Strip rotor