

Mini-6KT Mini Centrifuge



1. Mini-6KT Mini-Centrifuge can be supplied with different rotors for different applications.
2. It is also offered two-speed options (4,000 & 6,000rpm).
3. The electronic braking can reach quick deceleration to reduce handling time.

Feature:

1. Composite shell, shatter-resistant
2. The electronic braking system for rapid sample processing
3. Quick release rotor tub, allows for easy rotor exchange
4. Mini-6KT Two speed options: 4,000 & 6,000rpm
5. Patented rotor design for 48-well 50ul PCR plate Optional rotor: 4x5ml Eppendorf tube

Specification:

Model	Mini-6KT
Capacity	48well 50ul PCR rotor or 5 x 5ml tubes rotor
Two Speed options	4,000rpm & 6,000rpm
Maximum RCF	1,200 & 2,000 x g (6,000rpm)
Dimensions	135 x 147 x 120mm
Weight	0.7kg
Electrical	AC100~240V, 50/60Hz

Accessories:

Code	Description
AS-08100-00	Mini-6KT AC100V~240V 50~60Hz
AS-08101-01	48well 50ul PCR plate rotor
AS-08101-02	4 x 5ml tubes rotor (Optional)
AS-08101-03	Adaptors for 5ml normal tubes (Optional)