MiniT-3 Biological Indicator Incubator



Product Information

The MiniT-1 / MiniT-3 Biological Indicator Incubator are designed for temperature incubation of the EO Biological

Indicators. These provide easy verification of incubation temperature.

1. The Incubator have many wells for placement of the biological indicator vials and the wells can be custom-

made. User can set the incubating temperature.

2. Materials: Aluminum block

3. Well Diameter: According to custom-making

| Code | Description |
|-------------|-----------------------------------|
| AS-09020-00 | MiniT-1 AC100~240V, 50~60Hz, 0.5A |
| AS-09030-00 | MiniT-3 AC100~240V, 50~60Hz, 0.5A |

Features

1. LCD display

2. MiniT-1 is with DC12V power

MiniT-3 is with DC24V power

3. Operating Temperature Range:

MiniT-1: From 5°C above ambient to 80°C

MiniT-3: From 5°C above ambient to 100°C

3. Time Range: 1min ~ 99h59min

4. Temp. Control Accuracy: ≤±0.5°C

5. Display Accuracy: 0.1°C

6. Heating Time: < 6min (from 20°C to 57°C)

7. Temp. Control Accuracy @57°C: ≤±0.35°C

8. The Incubator have many wells for placement of the biological indicator vials and the wells can be custom-

made. User can set the incubating temperature and time.

9. Material: Aluminum block

Technical Data

| Model | MiniT-1 | MiniT-3 |
|------------------------------|---------------|---------------|
| Temperature Control Range | RT+5°C~80°C | RT+5°C~100°C |
| Timer Range | 1min~99h59min | 1min~99h59min |
| Temperature Control Accuracy | ±0.5°C | ±0.5°C |
| Display Accuracy | 0.1°C | 0.1°C |
| Heating Time (25°C to 57°C) | ≤6min | ≤6min |
| Temp. Control Accuracy @57°C | ±0.35°C | ±0.35°C |

| Dimensions (DxWxH) mm | 110 x 150 x 80 | 120 x 152 x 112 |
|-----------------------|----------------|-----------------|
| Net Weight | 0.5kg | 0.85kg |

Accessories

| Code | Туре | Description |
|-------------|-----------|-------------|
| AS-09021-01 | A-MiniT-1 | 15хФ10mm |
| AS-09021-02 | B-MiniT-1 | 12xФ12mm |
| AS-09021-03 | C-MiniT-1 | 6xФ15mm |
| AS-09021-04 | D-MiniT-1 | 6хФ16mm |
| AS-09021-05 | E-MiniT-1 | 24xФ8.8mm |
| AS-09021-06 | F-MiniT-1 | 18хФ9mm |