

# MiniT-H2C / MiniT-H3 Dry Bath Incubator



## Product Information

MiniT-H2C / MiniT-H3 are 3-in-1 mini dry bath incubator. The MiniT-H2C is a heating and cooling incubator and the MiniT-H3 is a heating only incubator with 3 individually controlled mini dry bath. Many aluminum blocks are available covering a range of tube sizes from 0.2ml, 0.5ml, 1.5ml or 2.0ml tubes blocks. These are suitable for labs with heavy work load.

## Features

1. Working area dimension: 3 individually controlled
2. LCD display. It is easy to setup and use
3. Simultaneously display the set and actual time, temperature, Auto-stop and sound-alert after program completed
4. MiniT-H3 has a heating function only and its temperature ranges from RT+5~100°C
5. MiniT-H2C has cooling/heating functions and its temperature ranges from 0°C~99.9°C
6. Temperature deflection adjusting meet with different demands

## Technical Data

Model	MiniT-H2C			MiniT-H3		
Working area	Block-A	Block-B	Block-C	Block-A	Block-B	Block-C
Temperature Control Range	0°C~99.9°C	RT.+5~99.9°C	RT.+5~99.9°C	RT.+5~99.9°C	RT.+5~99.9°C	RT.+5~99.9°C
Heating Time(25 to 99.9°C)	≤20min	≤20min	≤20min	≤20min	≤20min	≤20min
Cooling Time(25 to 0°C)	≤25min (*)	—	—	—	—	—
Temperature Control Accuracy	≤ 0.5°C			≤ 0.5°C		
Temp. Control Accuracy @37°C	0.2°C			0.2°C		
Display Accuracy	0.1°C			0.1°C		
Timing Range	1~999h or 1~999 min or 1~999sec			1~999h or 1~999 min or 1~999sec		

Block Quantity	3	3
Power	DC12V,120W	DC12V,120W
Dimension(mm)	360 × 180 × 150	360 × 180 × 150
Net Weight	3kg	3kg