# OS-200/OS-100



#### **Product Information**

The OS-100/OS-200 Orbital Shaker is a powerful variable speed shaker which provides efficient orbital motion. Speed and time are under microprocessor control.

The OS-200/OS-100 shaker can be applicable in different labs: microbiology, chemistry, immunology, biochemistry and molecular biology.

#### **Features**

- 1. Digital speed control
- 2. Speed range 50-250 RPM
- 3. Orbit 20 mm
- 4. Digital timer 1min 99h59min
- 5. Universal platform, flat platform or platforms with clamps (see below)
- 6. Maximum load capacity 2.5 kg
- 7. Orbital Shaker can be used in a cold room or incubator
- 8. Operating ambient: temperature range +5°C to +40°C
- 9. Overall dimensions: 280 x 270 x 110 mm
- 10.Weight: 7.2kg
- 11. Power supply: DC24V

### Platforms for OS-200

PP12-100 Platform with clamps for flasks, 100-150ml

PP06-250 Platform with clamps for flasks, 250-300ml

U P – 1 Universal platform with adjustable bars for different types of flasks, bottles

PP-1 Flat platform with non-slip rubber mat

### **Technical Data**

Model	OS-200	OS-100
Shaking Speed	50~250rpm	50~300rpm
Orbit	20mm	10mm
Independent timer	1min~99h59min	1min~99h59min
Load capacity	See below table	See below table
Load capacity	2.5kg	3.0kg
Motor	DC brushless	DC brushless
Operating ambient	4-45°C	4-45°C
Dimension(L x W x H) mm	280 x 270 x 110	280 x 270 x 110

#### **Accessories**

Interchangeable platforms for OS-100/OS-200

Code	Description
AS-04021-01	PP-1 Flat platform with non-slip rubber mat
AS-04021-02	PP12-100 Platform with clamps for flasks, 100-150ml
AS-04021-03	PP06-250 Platform with clamps for flasks, 250-300ml
AS-04021-04	UP-1 Universal platform with adjustable bars for different types of flasks, bottles
AS-04021-05	MP-1 Platform with springs

- 1. Digital speed control
- 2. Speed range 50-250 RPM
- 3. Orbit 20 mm
- 4. Digital timer 1min 99h59min
- 5. Universal platform, flat platform or platforms with clamps (see below)
- 6. Maximum load capacity 2.5 kg
- 7. Orbital Shaker can be used in a cold room or incubator
- 8. Operating ambient: temperature range +5°C to +40°C
- 9. Overall dimensions: 280 x 270 x 110 mm
- 10.Weight: 7.2kg
- 11. Power supply: DC24V

### Platforms for OS-200

- PP12-100 Platform with clamps for flasks, 100-150ml
- PP06-250 Platform with clamps for flasks, 250-300ml
- UP 1 Universal platform with adjustable bars for different types of flasks, bottles
- PP-1 Flat platform with non-slip rubber mat

## **Technical Data**

Model	OS-200	OS-100
Shaking Speed	50~250rpm	50~300rpm
Orbit	20mm	10mm
Independent timer	1min~99h59min	1min~99h59min
Load capacity	See below table	See below table
Load capacity	2.5kg	3.0kg
Motor	DC brushless	DC brushless
Operating ambient	4-45°C	4-45°C
Dimension(L x W x H) mm	280 x 270 x 110	280 x 270 x 110

## **Accessories**

Interchangeable platforms for OS-100/OS-200

Code	Description
AS-04021-01	PP-1 Flat platform with non-slip rubber mat
AS-04021-02	PP12-100 Platform with clamps for flasks, 100-150ml
AS-04021-03	PP06-250 Platform with clamps for flasks, 250-300ml
AS-04021-04	UP-1 Universal platform with adjustable bars for different types of flasks, bottles
AS-04021-05	MP-1 Platform with springs