

# TD2-12 ID-Centrifuge



## Product Information

The TD2-12 ID-Centrifuge holds automatic balance control and user can set speed and centrifugation time. It can be suitable many types ID-CARDS from different manufacturers.

1. Capacity: 12 Gel cards
2. Speed: Two stages (It can be variable according to user need)
3. Max. Speed: 1800rpm
4. Centrifugal Force: 54~272g
4. Timer: 1~20min
5. Power: AC100-240V, 50/60Hz, 85W
6. Dimensions: DxWxH=440x340x170 mm
7. Weight: 7.0kg

## Features

1. Silent centrifuge for up to 12 ID-cards.
2. Removable and exchangeable spinning head
3. Automatic balance control

4. Pre-defined centrifugal parameters for optimum Gel cards processing; acceleration, speed, braking and time
5. Digital interface visor control and visualization of centrifuge speed and time remaining for processing
6. Programmable audible alarms for end of centrifuge time periods

### Technical Data

Model	TD2-12
Capacity	12 Gel Cards
Speed Range	800~1800rpm
Centrifugal Force	54~272xg
Timer	1~20min
Noise	≤55db
Power	AC100~240V, 85W
Dimensions (DxWxH)	440 x 340 x 170mm
Net Weight	7kg