## TD2-24 ID Card Centrifuge



## **Product Information**

The TD2-24 ID-Centrifuge holds automatic balance control and user can set speed and centrifugation time. It can be suitable many types ID-CARDS from different manufacturers.

- 1. Capacity: 24 Gel cards
- 2. Speed: Two stages (It can be variable according to user need)
- 3. Max. Speed: 1500rpm
- 4. Centrifugal Force: 33~294g
- 4. Timer: 1~20min
- 5. Power: AC100-240V, 50/60Hz, 85W
- 6. Dimensions: DxWxH=530x410x180 mm
- 7. Weight: 11.5kg

## Features

- 1. Silent centrifuge for up to 24 ID-cards.
- 2. Removable and exchangeable spinning head
- 3. Automatic balance control

4. Pre-defined centrifugal parameters for optimum Gel cards processing; acceleration, speed, braking and time

5. Digital interface visor control and visualization of centrifuge speed and time remaining for processing

6. Programmable audible alarms for end of centrifuge time periods

## **Technical Data**

Model	TD2-24
Capacity	24 Gel Cards
Speed Range	500~1500 rpm
Centrifugal Force	33~294g
Timer	1~20min
Noise	≤55db
Power	AC100~240V, 85W
Dimensions (DxWxH)	530 x 410 x 180mm
Net Weight	11.5kg