

MS7-S

7" Square Magnetic Stirrer

Features

- Speed range of 0-1500rpm
- Max. stirring quantity of H₂O at 10L
- Ceramic work plate provides excellent chemical resistant performance



MS-S

Magnetic Stirrer

Features

- Maintenance free brushless DC motor
- Speed range of 0-1500 rpm
- Max. stirring quantity of H₂O at 20L
- Stainless steel work plate cover with ceramic material provides good chemical resistance performance



F100 protective cover

for non-hotplate model,
use with MS-S



Chemical resistance



Strong magnetism



Specifications

Work plate dimension[WxD]	
Work plate material	
Motor type	
Stirring positions	
Max. stirring quantity[H ₂ O]	
Max. magnetic bar[length]	
Speed range	
Speed display	
Protection class	
Motor input power	
Motor output power	
Power	
Voltage	
Dimension[WxDxH]	
Weight	
Permissible ambient temperature and humidity	

MS7-S

184x184mm(7 inch)

Glass ceramic

Shaded pole motor

1

10L

80mm

0-1500rpm

Scale

IP21

15W

1.5W

30W

100-120/200-240V,50/60Hz

215x360x112mm

3.8kg

5-40°C, 80%RH

MS-S

φ135mm(5 inch)

Stainless steel cover with ceramic

Brushless DC motor

1

20L

80mm

0-1500rpm

Scale

IP42

18W

10W

30W

100-240V,50/60Hz

160x280x85mm

2.8kg

5-40°C, 80%RH