

MX-T6-Pro LCD Digital Tube Roller

Features

- LCD displays the speed and time
- Speed range of 10-70rpm
- Six roller design for gentle rolling and easy to clean
- Continuous or timer operation



MX-T6-S Tube Roller

Features

- Speed range of 0-70rpm
- Six roller design for gentle rolling and easy to clean

DLAB tube rollers provide a gentle and highly efficient rolling motions, ideal for mixing blood samples, viscous substances and liquid-solid suspensions. Used in a variety of applications, such as the prevention of blood coagulation, immune precipitation, etc.



LCD display



Timer



Specifications

Amplitude
Motor type
Rocking motion
Max. load capacity
Number of rollers
Roller size[length]
Speed range
Speed display
Timer
Timer display
Time setting range
Operation mode
Voltage
Power
Dimension [DxWxH]
Weight
Permissible ambient temperature and humidity
Protection class

MX-T6-Pro	MX-T6-S
24mm	24mm
DC motor	DC motor
Rocking and rolling	Rocking and rolling
4kg	4kg
6	6
280mm	280mm
10 -70rpm	0-70rpm
LCD	Scale
Yes	-
LCD	-
1-1199min	-
Continuous/timer	Continuous
100-240V,50/60Hz	100-240V,50/60Hz
30W	25W
450x260x120mm	450x260x120mm
5.1kg	4.5kg
5-40°C,80%RH	5-40°C,80%RH
IP21	IP21