

SK-D3309-Pro LCD Digital 3D Shaker

Features

- Three-dimensional rocking motion with an angle of 9°
- Max. load capacity is 5kg with platform
- Speed range of 10-70rpm
- Double LCDs for independent display of time and speed
- Various of platform accessories



Timer



Wide range of accessories

SK-D1807-S LED Digital 3D Shaker

Features

- Three-dimensional rocking motion with the angle of 7°
- Max. load capacity is 3kg with platform
- Speed range of 10-80rpm
- Double LEDs for independent display of time and speed
- Various of platform accessories





Specifications

	SK-D3309-Pro	SK-D1807-S
Tilt angle	9°	7°
Max. load capacity (with platform)	5kg	3kg
Motor type	DC motor	DC motor
Motor input power	40W	16W
Motor output power	24W	10W
Speed range	10-70rpm	10-80rpm
Speed display	LCD	LED
Timer display	LCD	LED
Time settings range	1min-19h59min	1min-19h59min
Voltage	100-240V, 50/60Hz	100-240V, 50/60Hz
Power	50W	20W
Weight	9kg	2.7kg
Permissible ambient temperature and humidity	5-40°C, 80%RH	5-40°C, 80%RH
Dimension[WxDxH]	360x430x106mm	270x330x130mm
Protection class	IP21	IP21
Data connector	-	-



SK330.5



Cat. No.

Description

SK-O330-Pro
SK-L330-Pro
SK-R330-Pro
SK-D3309-Pro

SK-O180-Pro/S
SK-L180-Pro/S
SK-R1807-S
SK-D1807-S

18900155

Dish platform with non-slip mat and fixed string (8pcs)
Dimension[WxD]: 30x28.5cm



18202835

Silicone mattress for 16x 5~15ml tubes, use with SK330.5



18202846

Silicone mattress for 6x 50ml tubes, use with SK330.5



SK180.4



18900185

Dish attachment with non-slip mat
Dimension[WxD]: 24x21.5cm



SK180.5



18900200

Large size dish attachment with non-slip mat
Dimension[WxD]: 30x23cm



SK180.6



18900326

Dish attachment with non-slip mat
Dimension[WxD]: 28x28cm



18202852

Silicone mattress for 10x 5~10ml tubes



18202853

Silicone mattress for 4x 50ml tubes



SK330.1



18900027

Universal attachment with 4 bars
Size: 32x32cm
Capacity for flask: 25mlx21pcs,
50mlx15pcs, 100mlx12pcs,
250mlx6pcs, 500x4pcs



SK330.3



18900040

Lengthways roller attachment
Size: 32x32cm
Capacity for flask: 25mlx21pcs,
50mlx15pcs, 100mlx12pcs,
250mlx6pcs, 500x4pcs



SK180.1



18900025

Universal attachment with 3 bars
Size: 25x24.5cm
Capacity for flask: 25mlx10pcs,
50mlx8pcs, 100mlx6pcs,
250mlx4pcs, 500x2pcs



SK180.3



18900038

Lengthways roller attachment
Size: 25x24.5cm
Capacity for flask: 25mlx10pcs,
50mlx8pcs, 100mlx6pcs,
250mlx4pcs, 500x2pcs



18900036

Bar for SK330.1 universal attachment



18900037

Bar for SK180.1 universal attachment



18900280

Connect sticks for double dishes shaker





SK330.2



SK180.2



Cat. No.	Description	SK-O330-Pro SK-L330-Pro SK-R330-Pro SK-D3309-Pro	SK-O180-Pro/S SK-L180-Pro/S SK-R1807-S SK-D1807-S
18900028	Fixing clip attachment Size: 32 x 32cm	■	
18900026	Fixing clip attachment Size: 25 x 24cm		■
18900029	Fixing clip for flask volume 25mL, used with SK330.2[56pcs] SK180.2[36pcs]	■	■
18900030	Fixing clip for flask volume 50mL, used with SK330.2[28pcs] SK180.2[18pcs]	■	■
18900031	Fixing clip for flask volume 100mL, used with SK330.2[28pcs] SK180.2 [18pcs]	■	■
18900032	Fixing clip for flask volume 200/250 mL, used with SK330.2[13pcs] SK180.2[8pcs]	■	■
18900033	Fixing clip for flask volume 500mL, used with SK330.2[9pcs] SK180.2[8pcs]	■	■