

## dFlow

### Electronic Bottle-Top Dispenser

#### Features

- Digital operation
- Electronic control provides a fatigue-free operation
- Remote control panel prevents manual disturbance during operation
- Two operation modes:
- Traditional Dispensing, the liquid is continuously dispensed
- Stepper Function, for dispensing series



#### Specifications

<b>Volume Range</b>	0.1-99.9mL, Max piston lift is 10mL, resolution is 100 $\mu$ L
<b>Volume Accuracy</b>	Dispenser: R = 10mL $\pm$ 30 $\mu$ L; CV = 10ml $\pm$ 10 $\mu$ L Stepper: R = 1mL $\pm$ 6 $\mu$ L; CV = 1ML $\pm$ 9 $\mu$ L
<b>Operating Temperature Range</b>	10°C-30°C
<b>Quality Control</b>	DIN EN ISO 8655
<b>Contro Type</b>	External control