

Shaking Rotator

Operation method:	3D circling rotation	Input power:	6W
Angle of tilt:	20°	Working plate size:	240×170mm
Max. Loading weight:	1.5kg	Protective class:	IP21
Motor type:	AC motor	Permissible ambient temperature:	5-50°C
Speed range:	20rpm	Permissible relative moisture:	80%
Running mode:	Continuous	Dimensions:	200×160×160mm
Input voltage(frequency):	AC220V±10%(50/60Hz)	Net. Weight:	2.5Kg

>>> KJ-201A Microplate Shaker (Mini-Shaker)

Features:

- Elastic silicon plate, easy to clean and wear-resisting. It can fit various kinds of cell culture plates.
- High performance & low noise AC motor to make rapid and thorough mixing.
- Adjustable knob design can ensure stepless speed adjustment, and constant run within 60min.



>>> KJ-201B Microplate Shaker

Features:

- Open-type plate, easy to clean on the surface, 4 pcs of cell culture plates can be mixed at one time.
- High performance & low noise AC motor to make rapid and thorough mixing.
- Adjustable knob design can ensure stepless speed adjustment, and constant run within 60 min.

>>> KJ-201C Microplate Shaker

Features:

- Elastic silicon plate, easy to clean and wear-resisting. It can fit various kinds of cell plates. 6pcs of cell culture plates can be mixed at one time. High performance & low noise AC motor to make rapid and thorough mixing.
- Adjustable knob design can ensure stepless speed adjustment, and constant run within 60 min.



Technical Data:	KJ-201A	KJ-201B	KJ-201C
Operation method:	Circling	Circling	Circling
Circling diameter:	8.5mm	3.2mm	9mm
Rating speed:	800-2600rpm	600-2200rpm	200-2800rpm
Max. loading:	1 pcs microplate	2/4 pcs microplate	1-6 pcs microplate
Running mode:	Timing/Continuous	Timing/Continuous	Timing/Continuous
Timing range:	0-60min	0-60min	0-15min
Permissible relative moisture:	80%	80%	80%
Input voltage(frequency):	AC 220V±10% (50/60Hz)	AC 220V±10% (50/60Hz)	AC 220V±10% (50/60Hz)
Output power:	20W	10W	10W
Protective class:	IP21	IP21	IP21
Permissible ambient temperature:	5-50°C	5-50°C	5-50°C
Dimensions:	185×185×135mm	240×250×135mm	310×280×225mm
Net. Weight:	1.5Kg	3.8Kg	6.2Kg