


# Roller Mixer

## >>> KJMR-IV Rotator Mixer

### Features:

- 360 degree mild rotation mixing can keeps sample in suspension. Adapted to mixing of blood, urinary samples etc.
- Driven by High performance DC declaration Motor, Running stable, No noise, Long lifetime.
- Compact and reasonable structure makes it running stable. Adapted to various environment.
- Four types elastic clip could be randomly combined for use. Adapted to 1.5ml, 15ml, 50ml centrifuge tubes and test tubes with Ø8, Ø9, Ø12, Ø13, Ø16mm diameter.



## >>> KJMR-IVA Rotator Mixer


### Features:

- 360 degree mild rotation mixing can keeps sample in suspension. Adapted to mixing of blood, urinary samples etc.
- Driven by High performance DC declaration Motor, Running stable, No noise, Long lifetime.
- Compact and reasonable structure makes it running stable. Adapted to various environment.
- LED with backlight can display the current speed, time and operation mode, convenient for reading and setting.
- Precise digital time and speed setting can easily carry out the continuous operation mode.
- Wide-range voltage design. Avoid the interference to the equipment by instable voltage.
- Four types elastic clip could be randomly combined. Adapted to 1.5ml, 15ml, 50ml centrifuge tubes and test tubes with Ø8, Ø9, Ø12, Ø13, Ø16mm diameter.



Technical Data:	KJMR-IV	KJMR-IVA
Operation method:	360°reverse rotation	360°reverse rotation
Tilt angle of the plate:	38°	38°
Speed range:	6-32rpm	5-35rpm
Timing range:	—	0-999min
Clip Qty:	Large size 12pcs Medium size 24pcs Small size 24pcs	Large size 12pcs Medium size 24pcs Small size 24pcs
Running mode:	Continuous	Timing/Continuous
Plate diameter:	Ø300mm	Ø300mm
Input voltage(frequency):	AC 100-240V (50/60Hz)	AC 100-240V (50/60Hz)
Input power:	10W	10W
Protective class:	IP31	IP31
Permissible ambient temperature:	5-50°C	5-50°C
Permissible relative moisture:	80%	80%
Dimensions:	310×300×420mm	310×300×420mm
Net. Weight:	4.3Kg	4.0Kg

### KJMR-IV / KJMR-IVA Available plates:




Large size clip, for 50ml centrifuge tubes.



Medium size clip, for to Ø16 test tubes and 10ml, 15ml centrifuge tubes.



Small size clip: for Ø13 test tube, Ø12 test tube, 1.5ml centrifuge tube.



Extra small size clip: for Ø8-9mm test tube and 0.5ml centrifuge tube.