

2000系列多用途生物显微镜

THE 2000 SERIES MULTI-PURPOSE BIOLOGICAL MICROSCOPE



PH-2000D

PH-2000E

PH-2000SM

PH-2000V

技术参数 Specification	型号 Model				
	PH-2200D	PH-2000E	PH-2000SM	PH-2000V	
消色差物镜 Viewing head	4X、10X、40X(弹簧) 4X、10X、40X (spring)	•	•	•	•
目镜 Huygens eyepiecesng head	大视场目镜WF10X、(WF16X 选配) Wide field eyepieces WF10X、(WF16X optional)	•	•	•	•
观察头 Viewing head	单目头45°倾斜 Monocular,Inclined at 45°	•			
	铰链式双目头, 30°倾斜 Compensation Free Binocular Head Inclined at 30°		•		
	三目头 Trinocular head			•	
	示教头 Dual viewing head, one tube 30°, inclined and another tube vertical,360°rotatable				•
载物台 Stage	简易平台, 平台尺寸110mmX120mm Plain stage ,Stage size 110mmX120mm	•	•	•	•
聚光镜 Condenser	拨盘光栏 N.A.0.65 Condenser with disc diaphragm N.A.0.65	•	•	•	•
调焦 Focusing	粗微调焦 Coarse and fine Adjustment	•	•	•	•
电源 Power	220V 或 110V 交流输入 AC Input 220V or 110V,DC output 6V	•	•	•	•
照明系统 Illumination	LED 光源 (带充电选配) LED lamp (Chargeble optional)	•	•	•	•

60系列多用途生物显微镜

THE 60 SERIES MULTI-PURPOSE BIOLOGICAL MICROSCOPE



技术参数 Specification	型号 Model	
	XSP-60A	XSP-60B
消色差物镜 Viewing head	4X、10X、40X	•
放大倍率 Magnification factor	40X-1000X	•
目镜 Huygens eyepiecesng head	10X、25X广角目镜 10X、25Xwide-angle eyepiece	•
观察头 Viewing head	双目头 Binocular head	•
照明系统 Illumination	LED灯(上下光源) LED light (upper and lower light source)	•
电源供给 Power supply	3.0V (3节5号电池) 3.0V (3 no. 5 batteries)	•