

00系列单目生物显微镜

00 SERIES MONOCULAR BIOLOGICAL MICROSCOPE



XSP-03



XSP-04



XSP-05



XSP-06

技术参数 Specification		型号 Model			
		XSP-03	XSP-04	XSP-05	XSP-06
观察头 Viewing head	单目观察头 Monocular head	•	•	•	•
目镜 Huygens eyepiecesng head	H5X、H10X	•			
	H5X、H10X、H12.5X		•		
	H5X、H10X、H15X			•	
	H5X、H10X、H16X				•
消色差物镜 Achromatic objective	10X、40X、100X	•	•	•	•
载物台 Stage	单层平台，120mm×120mm，移动尺 Plain Stage with moving scale; Stage size:120mm*120mm	•	•	•	•
聚光镜 Condenser	阿贝式聚光镜，N.A=1.25带可变光栏和滤光片 N.A=1.25 Abbe condenser Iris diaphragm&fliter	•	•	•	•
粗调范围 Coarse adjustment range	50mm	•	•	•	•
微调范围 Fine adjustment range	1.8-2.2mm	•	•	•	•
照明 Illumination	平凹反光镜Φ50mm（外接电光源LED灯选配） Concave-plan reflector Φ50mm (extra led lam optional)	•	•	•	•

学生系列 显微镜

Students microscope series

性能 Microscope XSP-60

镜头采用优质光学玻璃，和专业医用显微镜一样，高清晰消色差镜

目镜：16倍 物镜：4倍 10倍 40倍

高低适中聚焦轮便于儿童操作

上下2个LED灯，室内室外通用，LED灯泡寿命长达5万小时

使用三节AA电池

16X Eyepiece 4X, 10X, 40X Objectives

40X, 100X & 400X Magnification

Real optical glass lenses

Dual focusing knobs

2 LED lights, above and below

Uses 3AA Batteries (not included)

