

100系列多用途生物显微镜

THE 100 SERIES MULTI-PURPOSE BIOLOGICAL MICROSCOPE



XSP-100D

XSP-100E

XSP-100SM

XSP-100V

技术参数 Specification	型号 Model			
	XSP-100D	XSP-100E	XSP-100SM	XSP-100V
消色差物镜 Viewing head	4X、10X、40X (弹簧)、100X (弹簧、油) 4X、10X、40X(S)、100X(S、Oil)			
目镜 Huygens eyepiecesng head	大视场目镜: WF10X、(WF16X选购件) Wide Field Eyepiece:WF10X、(WF16X optional)			
观察头 Viewing head	单目头 Monocular	•		
	铰链式双目头, 45°倾斜 Compensation Free Binocular Head Inclined at 45°		•	
	三目头 Trinocular head			•
	示教头 Dual viewing head, one tube 30°, inclined and another tube vertical, 360° rotatable			•
载物台 Stage	双层活动平台; 平台尺寸115X125mm Double Layer Mechanical Stage, Size 115X125mm			
调焦 Focusing	同轴粗微调, 调焦范围30mm, 微调格值0.02mm Coaxial Coarse and Fine Adjustment, Focusing Range 30mm, Focusing Interval 0.02mm			
聚光镜 Condenser	NA=1.25阿贝聚光镜带可变光栏和滤色片 Abbe NA=1.25 with Iris Diaphragm&Filter			
照明系统 Illumination	220V/6V20W或110V/6V20W卤素灯 (LED灯选配) Halogen Lamp 220V/6V20W or 110V/6V20W (LED lamp optional)			

104系列生物显微镜

THE 104 SERIES BIOLOGICAL MICROSCOPE



XSP-104

XSP-104V

技术参数 Specification	型号 Model	
	XSP-104	XSP-104V
消色差物镜 Viewing head	4X、10X、40X (弹簧)、100X (弹簧、油) 4X、10X、40X (S)、100X(S、Oil)	
目镜 Huygens eyepiecesng head	大视场目镜: WF10X、(WF16X选购件) Wide Field Eyepiece:WF10X、(WF16X optional)	
观察头 Viewing head	单目头 Monocular	•
	示教头 Dual viewing head, one tube 30°, inclined and another tube vertical, 360° rotatable	
载物台 Stage	双层活动平台; 平台尺寸115X125mm Double Layer Mechanical Stage, Size 115X125mm	
调焦 Focusing	同轴粗微调, 调焦范围30mm, 微调格值0.02mm Coaxial Coarse and Fine Adjustment, Focusing Range 30mm, Focusing Interval 0.02mm	
聚光镜 Condenser	NA=1.25阿贝聚光镜带可变光栏和滤色片 Abbe NA=1.25 with Iris Diaphragm&Filter	
照明系统 Illumination	白炽灯6V/20W, 220V或者110V (卤素灯6V/20W选配, LED灯选配) Incandescent lamp 6V/20W, 220V OR 110V (Halogen lamp 6V/20W optional, LED lamp optional)	