

学生系列显微镜

STUDENT SERIES MICROSCOPE

学生系列显微镜

STUDENT SERIES MICROSCOPE



XSP-60



XSP-60 (2)



XSP-60 (3)



XSP-102



XSP-102V

技术参数 Specification	型号 Model			
	XSP-60	XSP-60(2)	XSP-60(3)	
消色差物镜 Viewing head	4X、10X、40X(弹簧) 4X、10X、40X (spring)	•	•	•
目镜 Huygens eyepiecesng head	大视场目镜WF10X、(WF16X 选配) Wide field eyepieces WF10X、(WF16X optional)	•	•	•
观察头 Viewing head	单目头45°倾斜 Monocular,Inclined at 45°	•	•	•
载物台 Stage	简易平台, 平台尺寸90mmX90mm Plain stage ,Stage size 90mmX90mm	•	•	•
聚光镜 Condenser	拨盘光栏 N.A.0.65 Condenser with disc diaphragm N.A.0.65	•	•	•
调焦 Focusing	粗动调焦 Coarse focusing	•	•	•
电源 Power	220V 或 110V 交流输入 AC Input 220V or 110V,DC ouput 6V	•	•	•
照明系统 Illumination	LED 光源 (带充电选配) LED lamp (Chargeble optinal)	•	•	•

技术参数 Specification	型号 Model		
	XSP-102	XSP-102V	
消色差物镜 Viewing head	4X、10X、40X(弹簧) 4X、10X、40X (spring)	•	•
目镜 Huygens eyepiecesng head	大视场目镜WF10X、(WF16X 选配) Wide field eyepieces WF10X、(WF16X optional)	•	•
观察头 Viewing head	单目观察头 Monocular head	•	
	示教头 Dual viewing head, one tube 30°,inclined and another tube vertical,360°rotatable		•
转换器 Nosepiece	三孔 Triple	•	•
载物台 Stage	简易平台带台带压片 Plain Stage with paired clips	•	
	简易平台带移动尺 Plain Stage with moving scale		•
调焦 Focusing	同轴粗微调, 调焦范围 20mm Coaxial Coarse And Fine Adjustment,Focusing Range 20mm	•	•
聚光镜 Condenser	可变光栏 Variable light barrier	•	•
电源 Power	白炽灯6V/20W, 220V或者110V (卤素灯6V/20W 选配, LED 灯选配) Incandescent lamp 6V/20W,220V OR 110V (Halogen lamp 6V/20W optional,LED lamp optional)	•	•