

多用途生物显微镜

MULTI-PURPOSE BIOLOGICAL MICROSCOPE



XSP-113RT



XSP-115RT

技术参数 Specification	型号 Model	
	XSP-113RT	XSP-115RT
消色差物镜 Viewing head	4X、10X、40X	•
放大倍率 Magnification factor	40X-1000X	•
目镜 Huygens eyepiecesng head	大视场目镜: WF10X、25X广角目镜 Wide Field Eyepiece:WF10X、Wide-angle eyepiece	•
观察头 Viewing head	单目头 Monocular	•
	双目头 Binocular head	•
照明系统 Illumination	LED灯(上下光源) LED light (upper and lower light source)	•
电源供给 Power supply	3.0V (3节5号电池) 3.0V (3 no. 5 batteries)	•

视频显微镜

VIDEO MICROSCOPE



技术参数 Specification	型号 Model	
	视频显微镜 Video Microscope	
消色差物镜 Viewing head	4X、10X、40X(弹簧) 4X、10X、40X (spring)	•
目镜 Huygens eyepiecesng head	大视场目镜WF10X、(WF16X 选配) Wide field eyepieces WF10X、(WF16X optional)	•
观察头 Viewing head	显示屏 Display screen	•
载物台 Stage	简易平台带移动尺, 平台尺寸90mmX90mm Plain stage with moving scale, Stage size 90mmX90mm	•
聚光镜 Condenser	拨盘光栏 N.A.0.65 Condenser with disc diaphragm N.A.0.65	•
调焦 Focusing	粗微调焦 Coarse and fine Adjustment	•
电源 Power	220V 或 110V 交流输入 AC Input 220V or 110V, DC ouput 6V	•
照明系统 Illumination	LED 光源 (带充电选配) LED lamp (Chargeble optinal)	•

