

116系列 生物显微镜

The 116 series biological microscope



XSP-116L

116系列生物显微镜

THE 116 SERIES BIOLOGICAL MICROSCOPE



XSP-116B



XSP-116D



XSP-116L

	技术参数 Specification	型号 Model		
		XSP-116B	XSP-116D	XSP-116L
消色差物镜 Viewing head	4X、10X、40X(弹簧) 4X、10X、40X (spring)	•	•	•
目镜 Huygens eyepiecesng head	大视场目镜WF10X、(WF16X 选配) Wide field eyepieces WF10X、(WF16X optional)	•	•	•
观察头 Viewing head	铰链式双目头 Binocular head	•		
	单目头45°倾斜 Monocular, inclined at 45°		•	•
载物台 Stage	简易平台带移动尺, 平台尺寸90mmX90mm Plain stage with moving scale, Stage size 90mmX90mm	•	•	
	简易平台, 平台尺寸90mmX90mm Plain stage, Stage size 90mmX90mm			•
聚光镜 Condenser	拨盘光栏 N.A.0.65 Condenser with disc diaphragm N.A.0.65	•	•	•
调焦 Focusing	粗微调焦 Coarse and fine Adjustment	•	•	•
电源 Power	220V 或 110V 交流输入 AC Input 220V or 110V, DC output 6V	•	•	•
照明系统 Illumination	LED 光源 (带充电选配) LED lamp (Chargeble optional)	•	•	•

