

116系列生物显微镜

THE 116 SERIES BIOLOGICAL MICROSCOPE



XSP-116NL



XSP-116SM



XSP-116V

技术参数 Specification	型号 Model			
	XSP-116NL	XSP-116SM	XSP-116V	
消色差物镜 Viewing head	4X、10X、40X(弹簧) 4X、10X、40X (spring)	•	•	•
目镜 Huygens eyepiecesng head	大视场目镜WF10X、(WF16X 选配) Wide field eyepieces WF10X、(WF16X optional)	•	•	•
观察头 Viewing head	单目头45°倾斜 Monocular,Inclined at 45°	•		
	三目头 Trinocular head		•	
	示教头 Dual viewing head, one tube 30°,inclined and another tube vertical,360°rotatable			•
载物台 Stage	简易平台,平台尺寸90mmX90mm Plain stage ,Stage size 90mmX90mm	•		
	简易平台带移动尺,平台尺寸90mmX90mm Plain stage with moving scale,Stage size 90mmX90mm		•	•
聚光镜 Condenser	拨盘光栏 N.A.0.65 Condenser with disc diaphragm N.A.0.65	•	•	•
调焦 Focusing	粗动调焦 Coarse Adjustment	•		
	粗微调焦 Coarse and fine Adjustment		•	•
电源 Power	220V 或 110V 交流输入 AC Input 220V or 110V,DC output 6V	•	•	•
照明系统 Illumination	LED 光源 (带充电选配) LED lamp (Chargeble optional)	•	•	•

100系列 多用途生物显微镜

The 100 series multi-purpose biological microscope



XSP-100SM