

107系列多用途生物显微镜

THE 107 SERIES MULTI-PURPOSE BIOLOGICAL MICROSCOPE



XSZ-107



XSZ-107BN



XSZ-107D

技术参数 Specification	型号 Model			
	XSZ-107	XSZ-107BN	XSZ-107D	
消色差物镜 Viewing head	4X、10X、40X (弹簧)、100X (弹簧、油) 4X、10X、40X(S)、100X(S、Oil)	•	•	•
目镜 Huygens eyepiecesng head	大视场目镜: WF10X、(WF16X选购件) Wide Field Eyepiece:WF10X、(WF16X optional)	•	•	•
观察头 Viewing head	铰链式双目头 Compensation Free Binocular Head	•		
	平移式双目头, 45°倾斜 Sliding Binocular Head Inclined at 45°		•	
	单目头 Monocular			•
载物台 Stage	双层活动平台; 尺寸132X140mm, 移动范围75X45mm Double Layer Mechanical Stage Size 132X140mm, Moving Range 75X45mm	•	•	•
调焦 Focusing	同轴粗微调, 调焦范围30mm, 微调格值0.002mm Coaxial Coarse and Fine Adjustment, Focusing Range 30mm, Focusing Interval 0.002mm	•	•	•
聚光镜 Condenser	NA=1.25阿贝聚光镜带可变光栏和滤色片 Abbe NA=1.25 with Iris Diaphragm&Filter	•	•	•
照明系统 Illumination	220V/6V20W或110V/6V20W卤素灯 (LED灯选配) Halogen Lamp 220V/6V20W or 110V/6V20W (LED lamp optional)	•	•	•

107系列多用途生物显微镜

THE 107 SERIES MULTI-PURPOSE BIOLOGICAL MICROSCOPE



XSZ-107SM



XSZ-107V



XSZ-N107

技术参数 Specification	型号 Model			
	XSZ-107SM	XSZ-107V	XSZ-N107	
消色差物镜 Viewing head	4X、10X、40X (弹簧)、100X (弹簧、油) 4X、10X、40X(S)、100X(S、Oil)	•	•	•
目镜 Huygens eyepiecesng head	大视场目镜: WF10X、(WF16X选购件) Wide Field Eyepiece:WF10X、(WF16X optional)	•	•	•
观察头 Viewing head	铰链式三目头 Compensation Free Trinocular Head	•		
	示教头 Dual viewing head, one tube 30°, inclined and another tube vertical, 360° rotatable		•	
	平移式双目头, 45°倾斜 Sliding Binocular Head Inclined at 45°			•
载物台 Stage	双层活动平台; 尺寸132X140mm, 移动范围75X45mm Double Layer Mechanical Stage Size 132X140mm, Moving Range 75X45mm	•	•	•
调焦 Focusing	同轴粗微调, 调焦范围30mm, 微调格值0.002mm Coaxial Coarse and Fine Adjustment, Focusing Range 30mm, Focusing Interval 0.002mm	•	•	•
聚光镜 Condenser	NA=1.25阿贝聚光镜带可变光栏和滤色片 Abbe NA=1.25 with Iris Diaphragm&Filter	•	•	•
照明系统 Illumination	220V/6V20W或110V/6V20W卤素灯 (LED灯选配) Halogen Lamp 220V/6V20W or 110V/6V20W (LED lamp optional)	•	•	•