

瓶口加液器

Bottle Top Dispenser



数字视窗
Digital window



齿轮调节
Gear adjustment



耐酸耐碱
Acid and alkali resistant



安全阀门保护
Safety valve protection



型号丰富
Various models



技术参数

TECHNICAL DATA

型号 Model	容量范围 Capacity Range	增量 Increment	测试容量 Test Capacity	系统误差 Systemic Error	随机误差 Random Error	重复误差 Repeat Error
DA-2ML	0.4-2ml	0.05	1ml	10ul	10ul	2ul
DA-5ML	1-5ml	0.1	2.5ml	25ul	25ul	5ul
DA-10ML	2-10ml	0.25	5ml	50ul	50ul	10ul
DA-30ML	5-30ml	0.5	15ml	125ul	125ul	25ul
DA-60ML	10-60ml	1	30ml	250ul	250ul	50ul