

磁力搅拌器 Magnetic Stirrer



强磁性 防跳子
Strong magnetic



LCD数显恒温
LCD display



陶瓷涂层, 耐腐蚀
Anti-corrosion ceramic plate



安全回路, 过温保护
Safety circuit over temperature protection

技术参数

TECHNICAL DATA

型号 Model	HS-12	HS-17	HS-19	HSC-19T
工作盘尺寸 Working plate size	12*12cm	17*17cm	19*19cm	19*19cm
工作盘材料 Working plate material	铝盘 Aluminum plate	铝盘 Aluminum plate	防腐蚀陶瓷盘 Aluminum plate	防腐蚀陶瓷盘 Aluminum plate
转速 Speed	0-1600rpm	0-1600rpm	0-1600rpm	200-2000rpm
功能 Function	加热搅拌 Heating & stirring	加热搅拌 Heating & stirring	加热搅拌 Heating & stirring	加热搅拌 Heating & stirring
加热功率 Heating power	180W	500W	600W	600W
工作盘温度 Working plate temperature	室温-280°C Room temperature-280	室温-300°C Room temperature-300	室温-350°C Room temperature-350	室温-380°C Room temperature-350
搅拌量 Stirring capacity	1L	3L	5L	5L