

数显加热磁力搅拌器

Digital Display Hotplate Stirrer



强磁性，防跳子
Strong magnetic field



防腐台面
Anti-corrosion ceramic plate



PT1000，精准控温
Accurate temperature control



超温保护
Over temperature protection

技术参数

TECHNICAL DATA

型号 Model	HS5S	工作盘 Working plate	φ135mm 陶瓷防腐盘 Anti-corrosion ceramic plate
电机类型 Control mode	无刷电机 Brushless motor	工作盘温度 Plate temperature	260°C
最大搅拌量 Max stirring volume	2L (H ₂ O)	加热功率 Power	400W
转速范围 RPM	50-1500rpm	温度传感器 Temperature sensor	PT1000 / 外置PT1000 Inside and outside PT1000
磁子尺寸 Stir bar size	20-40mm	电压 Voltage	AC110/220V
旋转方向 Mixing direction	向右 Right	产品尺寸 Product size	155*220*700mm
外壳等级保护 Shell protection level	IP42	净重 Net weight	0.7KG