

红外灭菌器

Infrared Sterilizer



5 - 7 秒杀菌
Sterilization only takes 5-7 seconds



厌氧环境可使用
Can be used in an anaerobic environment



无明火，更安全
No open-flame, more safe



极易清理
Easy to clean

技术参数

TECHNICAL DATA

型号 Model	DS-900	DS-900S
功率 Power	150W	240W
中心区最高温度 Max. temperature inside the heating hole	825°C±50°C	
加热温度 Heating temperature	两档加热 400°C / 850°C Two temperatures available 400°C/850°C	
最大消毒外径 Max. diameter of sterilizer goods	φ14mm	φ35mm
加热区总长 Length of heating zone	150mm	100m
升温时间 Warm-up time	10min	