

低速离心机

Low Speed Centrifuge



4种程序，一键切换
Built-in 4 programs, one key operation



LCD显示
LCD display



电子门锁
Electronic door lock



操作简单
Easy to use



故障自动报警
Automatic fault alarm

技术参数

TECHNICAL DATA

型号 Model	LC-6S	LC-8S	LC-12S	LC-8P
最高转速 Max rotation speed	500-5000r/min	500-5000r/min	500-4000r/min	500-4000r/min
最大相对离心力 Max relative centrifuge force	3074xg	3074xg	2220xg	2220xg
最大容量 Max capacity	6*10/15ml(离心管) centrifuge tube	8*10/15ml(离心管) centrifuge tube	12*10ml(离心管) centrifuge tube	8*15ml(采血管) blood collection tube
	6*5/10ml(采血管) blood tube	8*5/10ml(采血管) blood tube	12*5/10ml(采血管) blood tube	内置4种程序 CGF:高浓缩生长因子纤维蛋白 离心制造 PRP:富血小板血浆离心制造 APRF:高级血浆富含纤维蛋白 IPRF:可注射血浆富含纤维蛋白
转速精度 Accuracy	± 30 r/min			
定时范围 Timing range	10s ~ 99min59s			
整机噪音 Noise	≤ 60dB(A)			
电源 Voltage	AC220V / 110V-50Hz 5A			
整机功率 Power	130W			
外形尺寸 Dimensions	375*320*240mm L*W*H			