

低速离心机

Low Speed Centrifuge



程序记忆
Program memory



多种转子可选
Multi rotors optional



温度保护&过震保护
Temperature protection &
overshock protection



双电子锁
Double electronic lock

技术参数

TECHNICAL DATA

型号 Model	LC500			
转速范围 Speed	1000 - 4000 rpm			
步进 Increment	100rpm			
转盘容量 Rotor capacity	15ml*12	10ml*24	50ml*6	100ml*4
	离心管/Centrifuge tube	采血管/Blood tube	离心管/Centrifuge tube	离心管/Centrifuge tube
最大相对离心力 Max relative centrifugal force	2075xg	2165xg	2165xg	2250xg
定时范围 Timing range	10sec - 60min			
温度显示精度 Temperature display accuracy	0.1°C			
开盖保护 Open cover protection	是 / Yes			
功率 Power	100W			