

# 水分分析仪

## Moisture Balance



多种保护，更安全  
Multiple protections, more secure



两种校准模式  
Two calibration modes



高效卤素加热灯  
High-efficiency halogen heating lamp



多设备打印，标准RS232端口  
Multi-device printing, standard RS232 interface

### 技术参数

#### TECHNICAL DATA

型号 Model	MBA1002	MBA1002a	MBA100b2	MBA1003
最大容量 Max capacity			100g	
测量精度 Measuring accuracy	0.01g	0.005g	0.002g	0.001g
实际刻度分度值 Actual scale dividing value			1mg	
检定分度值 Verification scale dividing value			10mg	
水分测试精度 Moisture test accuracy		±0.5%样品 / sample ≥2g		
水分含量决定可读性 Moisture content determines readability		0.02%~0.1%样品 / sample ≥2g		
温度调节允差 Temp control adjustment tolerance			±1°C	
加热温度范围设定 Heating temp set range	60~160°C (以1°C调整 / adjust based on 1°C)			
加热时间范围设定 Heating time set range	0min~99min(以1分钟调整 / adjust based on 1min)			
分析方法 Analysis method	手动、定时、自动 Manual / Timing / Automatic			
水分含量测定范围 Moisture content determination range	0%~100%			