

# 加热多联磁力搅拌器

## Multi Position Magnetic Stirrer With Heating



可设置同步或异步运行  
Can be set to run synchronously or asynchronously



LCD显示  
LCD display



304不锈钢盘面，耐腐蚀  
304 stainless steel plate, corrosion resistant



加热均匀  
Heat more uniformly



无刷电机，持久耐用  
Durable brushless motor

### 技术参数

#### TECHNICAL DATA

型号 Model	MHS4Pro	MHS6Pro	MHS8Pro
控制 Control	异步控制 (一键同步) Asynchronous control		
搅拌点位数 Number of stirring points	4	6	8
最大搅拌量 Max capacity	4*500ml	6*500ml	8*500ml
工作盘尺寸 Working plate size	4-(110*110mm)	6-(110*110mm)	8-(110*110mm)
工作盘材质 Plate material	不锈钢304 304 Stainless steel		
工作盘最高温度 Plate temperature	120°C		
定时 Timing	23h / 59min / 59sec		
适配转子 Bar size	Z30mm	Z30mm	Z30mm
功率 Power	240W	360W	360W
电压 Voltage	AC100-220V		