

微孔板离心机

Microplate Centrifuge



LCD显示
LCD display



微孔板转子，处理量大
Rotor for microplate,
more working capacity



无刷电机，动力强劲
Brushless motor with strong power



可定时调速
Speed adjustable and timing

技术参数

TECHNICAL DATA

型号 Model	MPC-96
显示 Display	LCD显示 LCD display
最大相对离心力 Max relative centrifugal force	600xg
最高转速 Max rotation speed	2300rpm±5%
转盘容量 Rotor capacity	细胞培养板 / cell cultrue plate * 4
计时/定时 Timing	59min59sec
输入功率 Input power	50W