

M系列微量可调移液器

M Series Micropipette



量程范围广
Wide range



分液精准
Dispense accurately



技术参数

TECHNICAL DATA

| 型号 Model | 容量范围 Capacity Range | 测试容量 Test Capacity | 吸头 Tip | 增量 Increment | 精度误差 Accuracy Error | 重复误差 Repeat Error |
|-------------|------------------------|-----------------------|-----------|-----------------|------------------------|----------------------|
| M2.5 | 0.1 - 2.5ul | 1.25 | 10ul | 0.05ul | ±3.00% | ±3.00% |
| M10 | 0.5 - 10ul | 5 | 10ul | 0.1ul | ±1.50% | ±0.80% |
| M20 | 2 - 20ul | 10 | 200ul | 0.5ul | ±1.20% | ±0.40% |
| M50 | 5 - 50ul | 25 | 200ul | 0.5ul | ±0.90% | ±0.30% |
| M100 | 10 - 100ul | 50 | 200ul | 1ul | ±1.00% | ±0.15% |
| M200 | 20 - 200ul | 100 | 200ul | 1ul | ±0.70% | ±0.15% |
| M1000 | 100 - 1000ul | 500 | 1000ul | 5ul | ±0.80% | ±0.20% |
| M5000 | 1000 - 5000ul | 2500 | 5000ul | 50ul | ±0.60% | ±0.15% |
| M10000 | 2 - 10ml | 5ml | 10ml | 0.1ml | ±1.20% | ±0.20% |