

多道移液器

Multichannel Pipette



数字视窗
Digital window



8 / 12道可选
8/12 channels optional



适用于96孔板
Suitable for 96-well plates



液头可360°旋转
The tip can rotate 360 degrees



技术参数

TECHNICAL DATA

型号 Model	规格 Type	容量范围 Capacity Range	吸头 Tip	增量 Increment	测试容量 Test Capacity	允许系统最大误差 Allow maximum error
M8 - 10	8道 8 channels	0.5 - 10ul	10ul	0.1ul	5	2.50% 0.125
M8 - 50	8道 8 channels	5 - 50ul	200ul	0.5ul	25	1.50% 0.375
M8 - 300	8道 8 channels	50 - 300ul	300ul	5ul	150	1.00% 1.500
M12 - 10	12道 12 channels	0.5 - 10ul	10ul	0.1ul	5	2.50% 0.125
M12 - 50	12道 12 channels	5 - 50ul	200ul	0.5ul	25	1.50% 0.375
M12 - 300	12道 12 channels	50 - 300ul	300ul	5ul	150	1.00% 1.500