

# 顶置搅拌机

## Overhead Stirrer



启动平缓  
Smooth start



搅拌速度可调  
Adjustable speed



直流永磁电机  
DC permanent magnet motor



搅拌桨位置可调  
Position of stirring paddle is adjustable

### 技术参数

#### TECHNICAL DATA

型号 Model	OS15	OS30
最大搅拌量 (H <sub>2</sub> O) Max amount of stirring(H <sub>2</sub> O)	10L	20L
转速范围 Range of rotate speed	100 - 2000rpm	50 - 1200rpm
定时 Timing	23h59min	23h59min
额定扭矩 Max torque	15Ncm	30Ncm
最大粘度 Max viscosity	10000mPas	10000mPas
转夹头直径范围 Clamping diameter range of drill chuck	1-8mm	1-8mm
额定功率 Power	40W	60W
电压 Voltage	AC 100-240V	AC 100-240V