

无油隔膜真空泵

Oilless Diaphragm Vacuum Pump



无油无污染
Oil-free and no pollution



带灯开关，安全美观
Illuminated switch



过载、过热、短路保护
Overload,overheat and short circuit protection



真空表精度高，反应快
High precision vacuum gauge

技术参数

TECHNICAL DATA

型号 Model	VP-10L	VP-15L
类型 Type	负压 Negative pressure	
抽气速率 Suction flow rate	10L/min	15L/min
极限真空度 Limit vacuum degree	0.075Mpa	0.08pa
噪音 Noise	<40db	<45db
外壳保护等级 Enclosure ratings	IP21	
电机功率 Power supply	20W	30W
隔膜阀片 Diaphragm	进口橡胶 Rubber	
电子保护 Protection devices	过载保护、过温保护、短路保护 Overload protection, over temperature protection, short circuit protection	
毛重/净重 Gross/net weight	1.95/1.58 KG	