

多管混匀仪

Multi Tube Vortex Mixer



软启动，均匀加速
Smooth start, stable operation



处理量大
Large handling capacity



无刷电机，动力强劲
Brushless motor with strong power



点动/连续，方便切换
Touch / Run, easy to switch

技术参数

TECHNICAL DATA

型号 Model	VM800
运行方式 Operating mode	点动 / 连续 Touch / Run
转速 Speed range	500 ~ 2800rpm
振幅 Amplitude	4mm水平回旋 4mm horizontal swing
时间设置 Time range	1min ~ 99h59min
最大承载量 Maximum load	4.5kg
输入功率 Power	100W

适配模块详情

ADAPTATION MODULE DETAILS

试管架型号 Model	试管架规格 Specs	适用规格 Model
VM801	50*φ10mm	φ9 ~ φ10mm
VM802 标配(Standard)	50*φ12mm	φ11 ~ φ13 mm (1.5/2ml离心管centrifuge tube)
VM803 标配(Standard)	50*φ15mm	φ14 ~ φ17mm (5/15ml离心管centrifuge tube)
VM804	50*φ20mm	φ18 ~ φ21mm
VM805	15*φ25mm	φ24 ~ φ27mm
VM806 标配(Standard)	15*φ29mm	φ28 ~ φ32 mm (50ml离心管centrifuge tube)